

TN-1396

Analysis of Ultrashort-chain PFAS in Tomatoes Using a Luna Polar Pesticides HPLC Column



Francesco Romaniello, PhD¹; Consolato Schiavone, PhD¹; Chiara Portesi, PhD¹; Luigi Margarucci, PhD²

¹Istituto Nazionale di Ricerca Metrologica, Strada delle Cacce 91, 10135 Turin

²Phenomenex s.r.l. Via Marino Serenzri 15D, 40013, Castel Maggiore Bologna

Overview

Food safety remains a top priority for both industry and consumers, with accurate contaminant quantification being essential to protect public health. Among these contaminants, per- and polyfluoroalkyl substances (PFAS) have garnered increasing attention due to their environmental persistence and links to serious health concerns, including cancer and immune system suppression. Short- and ultra-short-chain (USC) PFAS, in particular, pose unique challenges due to their high mobility and bioavailability in water-rich environments, making them more likely to enter the food chain especially through high-moisture foods.

While most current studies on USC-PFAS focus on simple matrices like drinking water, limited work has addressed their quantification in complex food samples. This study presents a targeted analytical method for detecting USC-PFAS (C1–C4, including TFA, PFPrA, PFBA, PFMeS, PFETs, PFPrS, PFBS, and DFA) in tomato for their high-water content and potential interferences from pigments and hydrophilic compounds. The newly developed chromatographic method improves retention and sensitivity for these analytes, offering a robust solution for PFAS detection in complex food matrices.

Sample Preparation

The extraction of tomato samples was conducted following the ion pair protocol described in EURL-POPs guidance [1].

Step	Description
1	Sample (2 g) + STD + ILIS in 50mL centrifuge tube
2	Add 2 mL Tetrabutylammoniumhydrogensulfate solution (0.5 M) + 4 mL Sodium carbonate solution (0.25 M) + 20 mL MTBE (Methyl tert butylether) + shaking (25 min)
3	Centrifugation (2500 x g, RT, 10 min)
4	Solvent evaporation to dryness at 40 °C
5	Dissolve residue in 300 µL MeOH/1% formic acid (2:1:v:v) in ultrasonic bath
6	Transfer supernatant in PP vial
7	UHPLC-MS/MS

LC Conditions

Column: Luna Polar Pesticides 3µm
Dimensions: 150 x 2.1mm
Delay Column: Luna Polar Pesticides 3µm
Dimensions: 100 x 2.1mm
Part No.: 00F-4798-AN
Mobile Phase: A (0.3% (v/v) HCOOH in H₂O)
 B (0.3% (v/v) HCOOH in ACN)
Gradient:

Time (min)	% B
0	10
1	10
2	50
13	90
14	100
15	100
16	10
22	10

Flow Rate: 0.400mL/min
Injection Volume: 10µL
Temperature: 45°C
LC System: SCIEX ExionLC™ AC System
Detection: Electrospray Ionization Tandem Mass Spectrometer (ESI-MS/MS)
Detector: QTRAP® 6500+

MS/MS Conditions

Ion Source: IonDrive® Turbo V
Polarity: Negative
Source Temperature: 400°C
GS1: 35psi
GS2: 45psi
CUR: 35psi
IS Voltage: -4000V

Table 1. MRM Transitions. *Quantifier ions

Analyte	Internal Standard	Q1 Mass (Da)	Q3 Mass (Da)	Dwell Time (msec)	DP (V)	EP (V)	CE (V)	CXP (V)
DFA_1*	13C_TFA	95.0	50.8	80	-50	-14	-19	-11
DFA_2	13C_TFA	95.0	95.0	80	-50	-14	-5	-11
TFA_1*	13C_TFA	113.0	68.9	80	-1	-3	-15	-6
TFA_2	13C_TFA	113.0	19.0	80	-1	-3	-40	-10
PFPrA_1*	13C3_PFPrA	163.0	118.9	80	-4	-14	-15	-10
PFPrA_2	13C3_PFPrA	163.0	19.0	80	-4	-14	-40	-10
PFBA_1*	13C3_PFBA	212.9	168.9	110	-23	-7	-14	-18
PFBA_2	13C3_PFBA	212.9	212.9	110	-23	-7	-5	-18
PFMeS_1*	13C4_PFBS	149.0	80.0	80	-20	-5	-31	-7
PFMeS_2	13C4_PFBS	149.0	98.9	80	-20	-5	-31	-7
PFETs_1*	13C4_PFBS	199.0	80.0	80	-42	-7	-32	-14
PFETs_2	13C4_PFBS	199.0	98.9	80	-42	-7	-32	-14
PFPrS_1*	13C4_PFBS	249.0	80.0	90	-17	-5	-65	-15
PFPrS_2	13C4_PFBS	249.0	98.9	90	-17	-5	-35	-15
PFBS_1*	13C4_PFBS	299.0	80.0	130	-50	-12	-73	-9
PFBS_2	13C4_PFBS	299.0	98.9	130	-50	-12	-40	-9

Figure 1. Chromatographic separation of ultrashort chain PFAS standards at 0.750 µg/kg for DFA, TFA, and PFPrA and 0.075 µg/kg for PFMeS, PFETs, PFPrS, PFBA, and PFBS




Table 2. USC-PFAS peak's identification and respective retention time.

Peak No.	Retention Time	Analytes
1	3.69	DFA
2	4.07	TFA
2	4.41	PFPrA
4	4.43	PFBA
5	4.88	PFMeS
6	4.92	PFETs
7	5.64	PFPrS
8	6.73	PFBS
9	8.95	PFETs from delay column
10	9.85	TFA from delay column

Figure 2. Chromatographic separation of TFA standard at 0.750 µg/kg (MRM transition 112.9/68.9).




Table 3. TFA peak's identification and respective retention time.

Peak No.	Retention Time	Analytes
1	4.07	TFA
2	9.85	TFA from delay column

Figure 3. Chromatographic separation of USC-PFAS standard in tomatoes matrix at 0.750 µg/kg for DFA, TFA, and PFPrA and 0.075 µg/kg for PFMeS, PFETs, PFPrS, PFBA, and PFBS




Table 4. Peak identification on the chromatogram

Peak No.	Analytes
1	DFA
2	TFA
2	PFPrA
4	PFBA
5	PFMeS
6	PFETs
7	PFPrS
8	PFBS

Calibration Curves

An isotope dilution method was employed to construct four distinct calibration curves, with seven concentration points each. Three of these curves were created in a solvent, while one involved spiking an extracted matrix. For DFA, TFA, and PFPrA, calibration standards were prepared at concentrations ranging from 0 to 1.500 µg/kg (specifically, 0, 0.015, 0.038, 0.075, 0.150, 0.375, 0.750, and 1.500 µg/kg). Each of these standards

contained a constant 0.375 µg/kg of isotopically labeled internal standards (specifically, ¹³C-TFA and ¹³C₃-PFPrA). Similarly, for PFMeS, PFETs, PFPrS, PFBA, and PFBS, calibration levels were established from 0 to 0.150 µg/kg (0, 0.002, 0.008, 0.015, 0.038, 0.075, 0.113, and 0.150 µg/kg), with an internal standard concentration of 0.075 µg/kg for ¹³C₃-PFBA and ¹³C₄-PFBS.

Figure 4. Calibration curves, average of five






Figure 4. Calibration curves, average of five






Matrix Effect and Method Precision

Method precision and Matrix effect: The method's precision was assessed by integrating the areas of labeled PFAS in calibration curves across 37 injections from different vials, both in solvent and matrix. The RSD% for ^{13}C -TFA, $^{13}\text{C}_3$ -

PFPrA, $^{13}\text{C}_3$ -PFBA, and $^{13}\text{C}_4$ -PFBS were 5.57, 2.74, 5.22, and 3.24, respectively (**Figure 5**).

Figure 5. Area of labeled USC-PFAS variation along 37 injection of standard in solvent and fortified matrix



Conclusion and Acknowledgment

This study presents a novel method for detecting ultra-short-chain PFAS in complex food matrices such as tomato products. It introduces a polar analytical and delay column setup that minimizes system interferences thanks to a huge delay of TFA contamination coming from system. This effect significantly improves sensitivity, reliability, and robustness.

In addition, the method shows excellent linearity ($R^2 > 0.99$) for all analytes, including TFA, in both solvent and matrix.

The authors acknowledge the invaluable support of the IMPreSA research infrastructure, funded through Regione Piemonte's INFRA-P POR-F.E.S.R 2014-2020 call for public research initiatives.

References

[1] C. Schiavone, F. Romaniello, P. Brizio, M. Suman, A. M. Rossi, C. Portesi, Reliable and sensitive ultra-short chain per- and polyfluorinated alkyl substances (PFAS) analysis in food: a polar reverse-phase chromatography approach, *Journal of Chromatography A*, 2025, 466332, ISSN 0021-9673, <https://doi.org/10.1016/j.chroma.2025.466332>

[2] European Union Reference Laboratory for POPs (EURL), Section 6.3.2, Module 1B: Solid-liquid extraction based on ion-pair extraction, in Annex V2.0 of the Guidance Document on Analytical Parameters for the Determination of Per- and Polyfluoroalkyl Substances (PFAS) in Food and Feed (issued 10 September 2024)



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

**BE-HAPPY™
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Phenomenex, the P icon, and the Phenomenex product and service marks mentioned herein are trademarks or registered trademarks of Phenomenex, Inc. in the United States and/or other countries. All other trademarks are the property of their respective owners. One or more authors are affiliated with Phenomenex, Inc. FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures. © 2025 Phenomenex, Inc. All rights reserved.

